

Guidescan2 Software Manual

Henri Schmidt

February 1, 2022

Contents

1	Installation	2
1.1	Software Dependencies	2
1.2	Build Procedure	2
2	Command Line Interface	3
2.1	Genome Indexing	3
2.2	Off-Target Enumeration	4
2.3	Kmer Set Generation	5
3	Off-Target Databases	5
3.1	Decoding	6
4	Guidescan2 Pipelines	6
4.1	Analysis of Individual gRNAs	6
4.2	Off-Target Database Construction	6
4.3	Off-Target Database Construction for F1-Cross Genomes	6
5	Miscellaneous Features	6

1 Installation

1.1 Software Dependencies

The off-target enumeration tool, Guidescan2, relies on the installation of a modern C++ compiler and associated buildtools. Explicitly, the tool has the following dependencies.

- CMake version 3.1.0 or higher
- C++ compiler that supports C++11 features such as:
 - regex
 - constexpr
 - default constructors
- C++ support for POSIX threads (pthreads)

It is known that GCC version 4.9 or greater and clang 3.1 or greater with libc++ are both supported. But to error on the side of caution, the most up to date and stable version of GCC is recommended. As often high-performance computing clusters do not possess recent versions of C++ compilers, it is therefore recommended that the software is built and executed in some sort of virtual environment. We have had success with Singularity and Docker for these special cases, though if a modern compiler is available, that can be used instead.

The aforementioned dependencies are the bare minimum to obtain anything useful out of the tool. However, it is strongly recommended, and necessary in some our pipelines, that biological computing tools are also available. In particular, the Samtools suite is useful since the output database will be in a SAM file format. There are also bindings for the tool in several programming languages that may help users to interface with the Guidescan2 databases programatically.

Helper scripts for various pipelines require a Python3 and Conda installation. We include a Conda environment file that contains all the necessary dependencies, though the scripts should run fine if standard data science packages are installed.

1.2 Build Procedure

Building the off-target enumeration tool, Guidescan2, is straightforward usign CMake once the dependencies have been installed. First set the working directory to the guidescan-cli project root. Then execute the following commands.

```
$ mkdir build; cd build
$ cmake -DCMAKE_BUILD_TYPE=RELEASE ..
$ make
```

The executable will be output in the `build/bin/` directory under the name `guidescan`. To make things convenient, we recommend that you add the program to your path. One way to do this is to place the executable under `$HOME/bin/` and then append this directory to your `$PATH`. The following code snippet does this and then permanently updates your `$PATH` in your `bashrc`.

```
$ mkdir -p $HOME/bin
$ cp bin/guidescan $HOME/bin
$ echo export PATH=$PATH:$HOME/bin >> ~/.bashrc
$ source ~/.bashrc
```

2 Command Line Interface

For ease of exposition, we will assume Guidescan2 is installed under your path and can be executed by simply running the command `guidescan`. However, everything will apply even if you execute Guidescan2 using its absolute path.

To run the Guidescan2 off-target enumeration and view its sub-commands, simply execute the program with the `--help` flag.

```
$ guidescan --help
```

There are three sub-commands for the Guidescan2 tool.

- **index**: Constructs a genomic index from a FASTA file.
- **enumerate**: Enumerates off-targets against a reference genomic index for a particular set of kmers.
- **http-server**: Spawns a local http-server that responds to off-target queries. That is, for a GET request containing a sequence of interest, enumerates all off-targets for that sequence.

The sub-commands each have their own interfaces, accessed as such:

```
$ guidescan [SUBCOMMAND] --help
```

For almost all use cases, it will suffice to stick to the main `enumerate` and `index` commands, the others are only helpers for very specific purposes.

2.1 Genome Indexing

```
$ guidescan index --help
```

The **index** sub-command takes a FASTA file as input and constructs a compressed genomic index that is used for off-target search. The index is stored as three files in the current working directory with the name ***.index.***. The index consists of a forward and reverse strand index as well as a small metadata file containing chromosome information. Auxiliary files with extensions **.forward.fna** and **.reverse.fna** are written to disk during the indexing process but these can be deleted once the genome is constructed.

The **index** step is the only step that requires a decent amount of memory to execute. For the human reference genome, 32GB of memory will suffice. However once the indices are constructed they can be transferred across devices as they are relatively small. Indices for a wide variety of genomes are available at <http://www.guidescan.org/>. And though this step requires a moderate amount of memory, it is quick to execute, taking under ten minutes on a standard laptop.

2.2 Off-Target Enumeration

```
$ guidescan enumerate --help
```

Once an index is constructed, off-target enumeration proceeds in an online fashion with the **enumerate** sub-command. Input to this command is a pointer to a reference index and a set of kmers to evaluate off-targets against. PAMs are included in this kmer set and are matched at the specified end of the kmer. The output is a database in SAM file format containing kmers that have passed the Guidescan2 filtering step. Complete off-target information is hex-encoded in the **of** field and can be decoded using an included Python3 script described later.

The command understands several options. Most importantly is the flag **-k/--mismatches** that describes the mismatch radius to search for kmer matches within. Typically, this parameter is set to a small value of, typically at most 3 or 4, as the search complexity grows exponentially with this parameter. Second, it is often useful to specify an alternative set of PAMs to match off-targets against since different CRISPR systems have varying PAM specificity. Of course this step can be emulated using repeated calls to **enumerate** with different PAMs, but we include this as a feature for convenience. Multi-threading across cores is available with the **--threads** option. This may result in out of order entries compared to the input kmer set since threads may complete off-target enumeration at different speeds.

We should also note that kmers that have multiple perfect matches are thrown away by default. This behavior can be turned off...

The **enumerate** command **requires** a set of genomic sequences to evaluate, we refer to these as the *kmer set*. The kmer set is specified in a very flexible CSV format described below. It can be automatically generated by the tool or manually specified. The kmer set consists of six columns:

```
id,sequence,pam,chromosome,position,sense
```

of which only **id**, **sequence** and **sense** are required, though all must be specified in the header. An example is shown in (Figure 1). Notice that for the third example, we don't specify a PAM. In this case, we simply match off-targets against the kmer with no regards to PAM specificity – off-targets don't need a PAM.

```
id,sequence,pam,chromosome,position,sense
ce11_example_1,GCCGTTTCAGGAGCTCGACGA,NGG,chrIV,5499,-
ce11_example_2,CAAAATATGAAATTTTCAAG,NGG,chrIV,23896,+
ce11_example_3,TCTACTGAAAGTTTGCAAAA,,chrIV,5499,-
ce11_example_4,TATAAACTGTCAAAGTTGAG,NGG,chrIV,23703,+
```

Figure 1: : An example kmer file for the *Caenorhabditis elegans* genome.

2.3 Kmer Set Generation

As mentioned earlier, *kmer set* generation can be performed manually, for example, when you want to evaluate a specific set of gRNAs. Additionally, Guidescan2 can also generate *kmer set* automatically. To generate *kmer sets* we include a Python3 script, **generate_kmers.py**, that generates a set of kmers against an input FASTA file. To use the script simply write

```
$ python generate_kmers.py FASTA_FILE
```

and the output CSV will be sent directly to **stdout**.

The script understands several options, though it is only mandatory to specify the FASTA file with which kmers are scanned. Additionally, one can specify PAM and kmer length, with defaults set to NGG and 20. By default kmers are only generated from chromosomes that are sufficiently long, but this can be specified using the **--min-chr-length** option. It also uses positional information as an identifier, but an additional prefix can be specified using the **--prefix** option.

The dependencies for the script are Python3 and the Biopython toolkit.

3 Off-Target Databases

Guidescan2 outputs databases in the SAM file format, a common file format used in Bioinformatics. At minimum, the Guidescan2 database contains the off-target set for each screened kmer as well as chromosomal structure information. Importantly, **the header is essential** for downstream processing; take care to not delete it accidentally. The additional fields position, chromosome, and strand are filled in only if they are specified in the *kmer set* for the corresponding kmer.

Though convenient for storage, compression, and a variety of bioinformatics tasks, the off-target information is difficult to pull out. Namely, off-target information is hex-encoded in the attributes

field with label `of` and is not human-readable. As such, we include a script to decode the off-target information into human-readable format.

3.1 Decoding

The script `decode_kmers.py` is a Python3 script that decodes Guidescan2 databases sequentially into a human-readable CSV format. It takes the Guidescan2 SAM/BAM database and the original FASTA file as input and outputs a CSV to `stdout`. The script outputs information in two modes, denoted *succinct* and *complete*, which we will describe below.

For most use cases *succinct* mode will be sufficient. In this mode each kmer’s off-targets are decoded and the information is summarized into multiple columns. Namely, we output the counts for matches at various distances and an overall score denoted specificity which is defined in the Supplementary Text. We note that the specificity is only defined for 20-mers with NGG PAMs; this column of the output will be empty otherwise. We also include additional metadata with obvious meaning upon examination.

In *complete* mode we explicitly write out the sequence and genomic location of all off-targets for each kmer. In addition, in the 20-mer case we include the CFD score from which our specificity score is derived. We warn that since some kmers can have up to hundreds of thousands of off-targets, *complete* mode can result in an extremely large output. If this is of concern we recommend that you re-run Guidescan2 with a smaller value of mismatch parameter and then decode the databases again.

To run the script simply write,

```
$ python scripts/decode_database.py GUIDESCAN_SAM_DB FASTA_FILE
```

The script **depends** on the two pickled files in the subdirectory `scripts/cfd/`. As such, it is required that the script and the sub-directory `scripts/cfd/` are in the location for the program to properly execute. By default this is the case, but we caution that when moving the script, move the `scripts/cfd/` sub-directory as well.

4 Guidescan2 Pipelines

4.1 Analysis of Individual gRNAs

4.2 Off-Target Database Construction

4.3 Off-Target Database Construction for F1-Cross Genomes

5 Miscellaneous Features